Sick Bay Bed
With 4 Drawers

The Sick Bay Bed with drawers gives school first aid rooms a comfortable place for students to rest and recover. The backrest adjusts to suit how a student needs to sit or lie down and four drawers underneath keep first aid supplies close at hand.

7 YRS



VE FURNITURE
DESIGNER SPACE SOLUTIONS

www.vef.com.au | 1300 966 640 | info@vef.com.au

Sick Bay Bed

Features & Highlights

- Three layers of foam cushion the bed and hold up to daily use.
- Adjustable backrest improves patient comfort.
- Medical grade vinyl upholstery resists common medical-grade disinfectants.
- Four, deep integrated drawers keep supplies organised.
- Optional locks are available for the drawers.
- E0 laminate base with steel legs forms a robust frame.

Technical Specifications

Model Number	Description	W	D	H	BH
SBB001	Left-hand side bedhead	1800	600	500	500
SBB001-MBL	Left-hand side bedhead- Mobile	1800	600	500	500
SBB002	Right-hand side bedhead	1800	600	500	500
SBB002-MBL	Right-hand side bedhead- Mobile	1800	600	500	500

Sick Bay Bed

Sick Bay Bed Range



Model: SBB001
Sick Bay Bed
Left-hand side bedhead
W1800 x D600 x H500
Bed Height: H500



Model: SBB001-MBL
Sick Bay Bed
Left-hand side bedhead
W1800 x D600 x H500
Bed Height: H500



Model: SBB002
Sick Bay Bed
Right-hand side bedhead
W1800 x D600 x H500
Bed Height: H500



Model: SBB002-MBL
Sick Bay Bed
Right-hand side bedhead
W1800 x D600 x H500
Bed Height: H500

E0 Laminate



Dark Oak S-004



Natural Oak S-005



Light Oak S-006



Grey S-008

White S-009

Luma Vinyl




Coral 2102




Deep Green 8477




Grass 8479



Ocean 8495



Baltic 8800



Terracotta 8471

Note: Materials are displayed for reference and may vary by screen.
For the truest look and feel, we recommend checking out a physical sample.

Price Group 3

Cleaning Guide

General Care Tips

Daily Care: Wipe the surface with a soft cloth dampened with warm water to remove dust, dirt and everyday marks. For routine cleaning, a mild detergent can be used to help lift grime while maintaining the finish of the product.

Deep Cleaning: For a more thorough clean, mix a small amount of mild dish soap with warm water and gently clean the surface using a soft cloth or sponge. Rinse with clean water to remove any residue and dry with a clean cloth to prevent streaks or water marks.

What to Avoid

Bleach and Harsh Chemicals: Avoid using bleach, acetone, solvent-based cleaners or harsh chemical products, as these may cause discoloration, surface damage or brittleness over time.

Abrasive Materials: Do not use steel wool, scouring pads or abrasive cleaning tools on polypropylene surfaces, as these can scratch and dull the finish.

VE FURNITURE
DESIGNER SPACE SOLUTIONS

www.vef.com.au | 1300 966 640 | info@vef.com.au